

Anti-P4Ha2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15460

Item Number: CM0026932



Introduction

P4HA2 Rabbit Polyclonal Antibody detects human P4HA2 (61kDa) via WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying collagen hydroxylation and extracellular matrix biology in heart, placenta, lung, and pancreas tissues.

[Learn More](#)

Parameter	Value
Product Name	P4HA2 Rabbit pAb
Remarks/Alias	MYP25; IncRNA-PE; P4HA2
Species	Human
Gene ID (Human)	8974
Gene ID	8974
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 239-333 of human P4HA2 (NP_004190.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15460
Protein Weight	61kDa
Shipping	Ice bag