

Anti-P4HA3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q7Z4N8

Item Number: CM0013194

Introduction

Rabbit polyclonal antibody to P4HA3 for WB, IF/ICC, ELISA. Detects human, mouse, rat P4HA3, a collagen prolyl 4-hydroxylase subunit. Suitable for extracellular matrix and fibrosis studies.

[Learn More](#)

Parameter	Value
Product Name	P4HA3 Rabbit pAb
Remarks/Alias	P4HA3
Species	Human
GeneID (Human)	283208
GeneID	283208
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-110 of human P4HA3 (NP_878907.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z4N8
Protein Weight	61kDa
Shipping	Ice bag