

Anti-P4Hb Rabbit Polymab® For Wb, Ip, If/Icc, Ihc-P, Elisa - P07237

Item Number: CM0002026



Introduction

P4HB Rabbit PolymAb® is a monoclonal antibody against protein disulfide-isomerase, validated for WB, IP, IF/ICC, IHC-P, and ELISA in human, mouse, and rat. Ideal for studying disulfide bond formation and ER protein folding.

[Learn More](#)

Parameter	Value
Product Name	P4HB Rabbit PolymAb®
Remarks/Alias	DSI; GIT; PDI; PHDB; PDIA1; PO4DB; PO4HB; PROHB; CLCRP1; ERBA2L; P4Hbeta
Species	Human
Gene ID (Human)	5034
Gene ID	5034
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 329-508 of human P4HB (NP_000909.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07237
Protein Weight	57 kDa
Shipping	Ice bag