

Anti-P53 Dinp1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96A56

Item Number: CM0012768

Introduction

Rabbit monoclonal antibody against human p53 DINP1, validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for autophagy and cancer research.

[Learn More](#)

Parameter	Value
Product Name	p53 DINP1 Rabbit mAb
Remarks/Alias	SIP; Teap; p53DINP1; TP53DINP1; TP53INP1A; TP53INP1B; p53 DINP1
Species	Human
Gene ID	94241
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human p53 DINP1 (Q96A56)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96A56
Protein Weight	27kDa
Shipping	Ice bag