


Anti-P53 (Phospho S392) Rabbit Monoclonal Antibody For Wb And Elisa - P04637

Item Number: CM0000757



Introduction

Phospho-specific rabbit monoclonal antibody against human p53 at serine 392, validated for Western blot and ELISA with rat cross-reactivity. Ideal for studying DNA damage response, apoptosis, and tumor suppression pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-S392 Rabbit mAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-S392
Species	Human
GeneID	7157
Immunogen	A synthetic phosphorylated peptide around S392 of human p53 (P04637)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag