


Anti-P53 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P04637

Item Number: CM0002946



Introduction

Rabbit monoclonal antibody against human p53 (P04637). KO-validated for Western blot, immunoprecipitation, and ELISA. Reacts with human. Ideal for tumor suppressor and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] p53 Rabbit mAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; 53
Species	Human
Gene ID (Human)	7157
Gene ID	7157
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-393 of human p53 (NP_000537.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag