


Anti-P53 Rabbit Monoclonal Antibody (Ko Validated) For Wb, Ihc-P, If/Icc - P02340

Item Number: CM0010622



Introduction

KO validated p53 rabbit monoclonal antibody for Western blot, IHC-P, and IF/ICC. Detects human, mouse, rat p53. Synthetic peptide immunogen, 43 kDa. For cancer research and apoptosis studies.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] p53 Rabbit mAb
Remarks/Alias	bbl; bfy; bhy; p44; p53; Tp53
Species	Mouse
GeneID (Human)	7157
GeneID	22059
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02340
Protein Weight	43 kDa
Shipping	Ice bag