


Anti-P53 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P02340

Item Number: CM0018896



Introduction

High-quality rabbit polyclonal antibody targeting p53 tumor suppressor. Validated for WB, IHC-P, IF/ICC, IP, and ELISA applications in human and mouse samples. Suitable for apoptosis, cell cycle, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	p53 Rabbit pAb
Remarks/Alias	bbl; bfy; bhy; p44; p53; Tp53
Species	Mouse
GeneID (Human)	7157
GeneID	22059
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 290-390 of mouse p53 (NP_035770.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02340
Protein Weight	43kDa
Shipping	Ice bag