


Anti-P53 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P04637

Item Number: CM0029942



Introduction

Rabbit polyclonal antibody against human p53, validated for Western blotting, immunoprecipitation, and ELISA. Detects the 44kDa p53 tumor suppressor protein, a key regulator of cell cycle arrest and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	p53 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMF55; TRP53; p53
Species	Human
Gene ID	7157
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 218-390 of human p53 (NP_035770.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag