

Anti-P53 Rabbit Polymab® (Ko Validated) For Wb, Ihc-P, If/Icc, Elisa - P04637

Item Number: CM0006844



Introduction

KO-validated rabbit monoclonal antibody against human p53. Suitable for WB, IHC-P, IF/ICC, and ELISA. Recombinant protein immunogen, cross-reacts with human. Ideal for cancer research and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] p53 Rabbit PolymAb®
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53
Species	Human
GenelD (Human)	7157
GenelD	7157
Immunogen	Recombinant protein of human p53.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag