

Anti-P60 Caf1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13112

Item Number: CM0010482



Introduction

Rabbit monoclonal antibody targeting human p60 CAF1 (CHAF1B), a chromatin assembly factor involved in nucleosome assembly during DNA replication and repair. Validated for WB, IHC-P, and ELISA applications.

[Learn More](#)

Parameter	Value
Product Name	p60 CAF1 Rabbit mAb
Remarks/Alias	CAF1; MPP7; CAF-1; CAF1A; CAF1P60; CAF-IP60; MPHOSPH7; p60 CAF1
Species	Human
Gene ID (Human)	8208
Gene ID	8208
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 459-559 of human p60 CAF1 (Q13112)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q13112
Protein Weight	61kDa
Shipping	Ice bag