

Anti-P60 Caf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13112

Item Number: CM0024710



Introduction

Rabbit polyclonal antibody against human p60 CAF1 (CHAF1B, CAF-1 p60) validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat. Ideal for chromatin assembly and DNA replication studies.

[Learn More](#)

Parameter	Value
Product Name	p60 CAF1 Rabbit pAb
Remarks/Alias	CAF1; MPP7; CAF-1; CAF1A; CAF1P60; CAF-IP60; MPHOSPH7; p60 CAF1
Species	Human
Gene ID (Human)	8208
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 270-559 of human p60 CAF1 (NP_005432.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13112
Protein Weight	61kDa
Shipping	Ice bag