

Anti-P63 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H3D4

Item Number: CM0017601



Introduction

Rabbit polyclonal antibody against human p63 (TP63) validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Research use for apoptosis, epithelial stem cell biology, and limb development.

[Learn More](#)

Parameter	Value
Product Name	p63 Rabbit pAb
Remarks/Alias	AIS; KET; LMS; NBP; RHS; p40; p51; p63; EEC3; OFC8; p73H; p73L; SHFM4; TP53L; TP73L; p53CP; TP53CP; B(p51A); B(p51B)
Species	Human
Gene ID (Human)	8626
Gene ID	8626
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human p63 (NP_003713.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H3D4
Protein Weight	45-77 kDa
Shipping	Ice bag