

## Anti-P63 Rabbit Monoclonal Antibody For Ihc-P, Elisa - Q9H3D4

Item Number: CM0011853



### Introduction

p63 rabbit monoclonal antibody for IHC-P and ELISA. Reacts with human, mouse, rat. Suitable for epithelial morphogenesis and stem cell biology research. SWISS: Q9H3D4.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | p63 Rabbit mAb                                                                                                           |
| Remarks/Alias    | AIS; KET; LMS; NBP; RHS; p40; p51; p63; EEC3; OFC8; p73H; p73L; SHFM4; TP53L; TP73L; p53CP; TP53CP; B(p51A); B(p51B)     |
| Species          | Human                                                                                                                    |
| Gene ID          | 8626 (Human)                                                                                                             |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 581-680 of human p63 (NP_003713.3). |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | IHC-P, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | Q9H3D4                                                                                                                   |
| Protein Weight   | 45-77 kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                                  |