

Anti-P63 Rabbit Monoclonal Antibody For Wb, Elisa, Chip - Q9H3D4

Item Number: CM0005135



Introduction

High-quality rabbit monoclonal antibody targeting human p63 (TP63). Validated for WB, ELISA, and ChIP, it detects human and rat p63 with a molecular weight of 45-77 kDa. Ideal for studying transcription regulation, epithelial morphogenesis, and limb development.

[Learn More](#)

Parameter	Value
Product Name	p63 Rabbit mAb
Remarks/Alias	AIS; KET; LMS; NBP; RHS; p40; p51; p63; EEC3; OFC8; p73H; p73L; SHFM4; TP53L; TP73L; p53CP; TP53CP; B(p51A); B(p51B)
Species	Human
GeneID (Human)	8626
GeneID	8626
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 581-680 of human p63 (NP_003713.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Rat
SWISS	Q9H3D4
Protein Weight	45-77 kDa
Shipping	Ice bag