

Anti-P63 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9H3D4

Item Number: CM0003385



Introduction

Recombinant rabbit monoclonal antibody targeting human p63 (Tumor protein 63). Validated for WB, IHC-P, and ELISA. Cross-reacts with human and mouse. Suitable for epithelial development and cancer research.

[Learn More](#)

Parameter	Value
Product Name	p63 Rabbit mAb
Remarks/Alias	AIS; KET; LMS; NBP; RHS; p40; p51; p63; EEC3; OFC8; p73H; p73L; SHFM4; TP53L; TP73L; p53CP; TP53CP; B(p51A); B(p51B)
Species	Human
Gene ID (Human)	8626
Gene ID	8626
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-203 of human p63 (Q9H3D4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H3D4
Protein Weight	45-77 kDa
Shipping	Ice bag