

Anti-P63 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9H3D4

Item Number: CM0019873



Introduction

Rabbit polyclonal antibody against human p63 (TP63), suitable for WB, IF-P, IHC-P, ELISA. Recognizes human, mouse, rat. Ideal for transcription regulation and epithelial morphogenesis research.

[Learn More](#)

Parameter	Specification
Product Name	p63 Rabbit pAb
Remarks/Alias	AIS, KET, LMS, NBP, RHS, p40, p51, p63, EEC3, OFC8, p73H, p73L, SHFM4, TP53L, TP73L, p53CP, TP53CP, B(p51A), B(p51B)
Species	Human
GeneID (Human)	8626
GeneID	8626
Immunogen	Synthetic peptide, A synthetic peptide corresponding to a sequence within amino acids 166-362 of human p63 (NP_001108452.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H3D4
Protein Weight	45-77 kDa
Shipping	Ice bag