


Anti-P70 S6 Kinase 1 Monoclonal Antibody For Wb, If/Icc, Elisa - P23443

Item Number: CM0011979



Introduction

High-quality Anti-p70 S6 Kinase 1 Rabbit monoclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets key mTOR downstream kinase regulating translation and cell growth. For IVD research.

[Learn More](#)

Parameter	Value
Product Name	p70 S6 Kinase 1 Rabbit mAb
Remarks/Alias	S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha; p70 S6 Kinase 1
Species	Human
GeneID	6198
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human p70 S6 Kinase 1 (NP_003152.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23443
Protein Weight	59 kDa
Shipping	Ice bag