

Anti-P70 S6 Kinase 1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P23443

Item Number: CM0021499



Introduction

Rabbit polyclonal antibody targeting human p70 S6 Kinase 1 (RPS6KB1). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse, rat. Ideal for mTOR signaling and cell proliferation research.

[Learn More](#)

Parameter	Value
Product Name	p70 S6 Kinase 1 Rabbit pAb
Remarks/Alias	S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha; p70 S6 Kinase 1
Species	Human
GeneID (Human)	6198
GeneID	6198
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 421-520 of human p70 S6 Kinase 1 (NP_003152.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23443
Protein Weight	59kDa
Shipping	Ice bag