

Anti-P70 S6 Kinase 2 Polyclonal Antibody For Wb, Elisa - Q9Ubs0

Item Number: CM0014524

Introduction

Rabbit polyclonal antibody targeting human p70 S6 Kinase 2 (RPS6KB2) for Western blot and ELISA applications. Cross-reacts with mouse and rat. Suitable for mTOR signaling, cell growth, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	p70 S6 Kinase 2 Rabbit pAb
Remarks/Alias	KLS; SRK; S6K2; S6KB; STK14B; S6KI(2); S6Kbeta; p70S6Kb; P70-beta; S6K-beta2; P70-beta-1; P70-beta-2; p70(S6K)-beta; p70 S6 Kinase 2
Species	Human
GeneID (Human)	6199
GeneID	6199
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human p70 S6 Kinase 2 (NP_003943.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBS0
Protein Weight	53kDa
Shipping	Ice bag