

Anti-P73 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O15350

Item Number: CM0006984



Introduction

Rabbit monoclonal antibody targeting human p73 (O15350). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat p73. p73 is a p53-related transcription factor involved in apoptosis, tumor suppression, and ciliogenesis.

[Learn More](#)

Parameter	Value
Product Name	p73 Rabbit mAb
Remarks/Alias	P73; CILD47; p73
Species	Human
Gene ID	7161
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human p73 (O15350).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15350
Protein Weight	70kDa
Shipping	Ice bag