

Anti-P75Ntr/Ngfr (Cd271) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P08138

Item Number: CM0003820



Introduction

High-specificity rabbit monoclonal antibody targeting human p75NTR/NGFR (CD271/TNFRSF16). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for neurotrophin signaling and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	p75NTR/NGFR Rabbit mAb
Remarks/Alias	CD271; p75NTR; TNFRSF16; p75(NTR); Gp80-LNGFR; p75NTR/NGFR
Species	Human
GeneID (Human)	4804
GeneID	4804
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 300-400 of human p75 NGF Receptor (P08138).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08138
Shipping	Ice bag