


Anti-P75Ntr/Ngfr Rabbit Polyclonal Antibody For Wb, Elisa - P08138

Item Number: CM0020352



Introduction

Rabbit polyclonal antibody against human p75NTR/NGFR, validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for neurotrophin signaling and neuronal differentiation research.

[Learn More](#)

Parameter	Value
Product Name	p75NTR/NGFR Rabbit pAb
Remarks/Alias	CD271; p75NTR; TNFRSF16; p75(NTR); Gp80-LNGFR; p75NTR/NGFR
Species	Human
GeneID (Human)	4804
GeneID	4804
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08138
Protein Weight	45kDa
Shipping	Ice bag