


Anti-P75Ntr/Ngfr Rabbit Polymab® For Wb, Ihc-P, Elisa - P08138 / Q9Z0W1

Item Number: CM0012125



Introduction

Rabbit monoclonal antibody targeting p75NTR/NGFR for western blot, IHC-P, and ELISA. Recognizes human and mouse NGFR with cross-reactivity to rat. Ideal for neurotrophin signaling, cell fate, and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	p75NTR/NGFR Rabbit PolymAb®
Remarks/Alias	CD271; p75NTR; TNFRSF16; p75(NTR); Gp80-LNGFR
Species	Human/Mouse
GeneID (Human)	4804
GeneID	4804/18053
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of human p75NTR/NGFR (NP_150086.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08138/Q9Z0W1
Protein Weight	45kDa
Shipping	Ice bag