

Anti-Pabpc1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P11940

Item Number: CM0015808



Introduction

Rabbit polyclonal antibody targeting human PABPC1, validated for western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for mRNA metabolism and translation research.

[Learn More](#)

Parameter	Value
Product Name	pAbPC1 Rabbit pAb
Remarks/Alias	PAB1; PABP; PABP1; PABPC2; PABPL1; pAbPC1
Species	Human
Gene ID (human)	26986
Gene ID	26986
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11940
Protein Weight	71kDa
Shipping	Ice bag