

Anti-Pabpc4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13310

Item Number: CM0028445

Introduction

Rabbit polyclonal antibody against human PABPC4. Applications: WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying mRNA poly(A) binding and post-transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	PABPC4 Rabbit pAb
Remarks/Alias	APP1, APP-1, PABP4, iPABP, PABPC4
Species	Human
GeneID (Human)	8761
GeneID	8761
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human PABPC4 (NP_003810.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13310
Protein Weight	71kDa
Shipping	Ice bag