

Anti-Pabpc5 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Du9

Item Number: CM0027003

Introduction

Rabbit polyclonal antibody targeting PABPC5 (polyadenylate-binding protein 5) for Western blot and ELISA applications. Demonstrates cross-reactivity with human and mouse species. Ideal for IVD research into mRNA metabolism and cytoplasmic regulatory processes.

[Learn More](#)

Parameter	Value
Product Name	PABPC5 Rabbit pAb
Remarks/Alias	PABP5; PABPC5
Species	Human
GeneID (Human)	140886
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-382 of human PABPC5 (NP_543022.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96DU9
Protein Weight	43kDa
Shipping	Ice bag