

Anti-Pabpn1 Polyclonal Antibody For Wb And Elisa - Q86U42

Item Number: CM0028652



Introduction

Human PABPN1 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects Human and Mouse PABPN1, a nuclear poly(A)-binding protein involved in mRNA 3'-end processing. Ideal for research in mRNA metabolism and muscular dystrophy.

[Learn More](#)

Parameter	Value
Product Name	PABPN1 Rabbit pAb
Remarks/Alias	OPMD; PAB2; PABII; PABP2; PABP-2; PABPN1
Species	Human
Gene ID	8106
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PABPN1 (NP_004634.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86U42
Protein Weight	33kDa
Shipping	Ice bag