

## Anti-Pabpn1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q86U42

Item Number: CM0003130



### Introduction

Reacts with human, mouse, rat PABPN1; suitable for Western blot, immunofluorescence, IHC, and ELISA; targets a nuclear poly(A)-binding protein involved in mRNA processing and muscle gene regulation.

[Learn More](#)

| Parameter        | Value                                                                                              |
|------------------|----------------------------------------------------------------------------------------------------|
| Product Name     | PABPN1 Rabbit mAb                                                                                  |
| Remarks/Alias    | OPMD, PAB2, PABII, PABP2, PABP-2, PABPN1                                                           |
| Species          | Human                                                                                              |
| Gene ID (Human)  | 8106                                                                                               |
| Gene ID          | 8106                                                                                               |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PABPN1 (Q86U42). |
| Source           | Rabbit                                                                                             |
| Category         | Monoclonal Antibodies                                                                              |
| Application      | WB, IF/ICC, IF-P, IHC-P, ELISA                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                  |
| SWISS            | Q86U42                                                                                             |
| Protein Weight   | 31kDa/33kDa                                                                                        |
| Shipping         | Ice bag                                                                                            |