


Anti-Pacsin2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Unf0

Item Number: CM0022742



Introduction

Rabbit polyclonal antibody against PACSIN2 (Syndapin-2), validated for WB and ELISA in human and mouse samples. PACSIN2 regulates caveolae morphogenesis, endocytosis, and HIV-1 spread.

[Learn More](#)

Parameter	Value
Product Name	PACSIN2 Rabbit pAb
Remarks/Alias	SDP11; PACSIN2
Species	Human
GeneID (Human)	11252
GeneID	11252
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-410 of human PACSIN2 (NP_009160.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UNF0
Protein Weight	56kDa
Shipping	Ice bag