

Anti-Pacsin3 Polyclonal Antibody For Wb And Elisa - Q9Uks6

Item Number: CM0030020

Introduction

Rabbit polyclonal antibody to PACSIN3, suitable for Western blotting and ELISA. Detects human, mouse, and rat PACSIN3, a protein involved in endocytosis and membrane protein trafficking.

[Learn More](#)

Parameter	Value
Product Name	PACSIN3 Rabbit pAb
Remarks/Alias	SDP111; PACSIN3
Species	Human
GeneID (Human)	29763
GeneID	29763
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 175-424 of human PACSIN3 (NP_057307.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKS6
Protein Weight	48kDa
Shipping	Ice bag