

Anti-Padi1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9ULc6

Item Number: CM0016368



Introduction

Rabbit polyclonal antibody against human PADI1, suitable for WB and ELISA. Cross-reacts with mouse and rat. Used for studying protein arginine deiminase in epidermal keratinocytes. Research use only.

[Learn More](#)

Parameter	Value
Product Name	PADI1 Rabbit pAb
Remarks/Alias	PDI; PADI1; PDI1; HPAD10; PADI1
Species	Human
GeneID (Human)	29943
GeneID	29943
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-277 of human PADI1 (NP_037490.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9ULC6
Protein Weight	75kDa
Shipping	Ice bag