

Anti-Padi1/Pad1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9ULc6

Item Number: CM0007636



Introduction

High-quality rabbit monoclonal antibody against human PADI1/PAD1, validated for WB and ELISA, cross-reactive with mouse. Ideal for studying protein citrullination and epidermal biology. Predicted molecular weight of 74kDa.

[Learn More](#)

Parameter	Value
Product Name	PADI1/PAD1 Rabbit mAb
Remarks/Alias	PADI1; HPAD10; PAD1; PDI; PDI1; peptidyl arginine deiminase 1; PADI1/PAD1
Species	Human
Gene ID (Human)	29943
Gene ID	29943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-277 of human PADI1/PAD1 (NP_037490.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9ULC6
Protein Weight	74kDa
Shipping	Ice bag