

Anti-Padi4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Um07

Item Number: CM0020523



Introduction

Rabbit polyclonal antibody targeting PADI4 (protein-arginine deiminase type-4). Verified for WB and ELISA. Reacts with human and mouse. Useful for histone code, chromatin regulation, and NETs studies.

[Learn More](#)

Parameter	Value
Product Name	PADI4 Rabbit pAb
Remarks/Alias	PAD; PAD4; PDI4; PDI5; PADI5; PADI4
Species	Human
Gene ID (Human)	23569
Gene ID	23569
Immunogen	synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 350-450 of human PADI4 NP_036519.2 .
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UM07
Protein Weight	74kDa
Shipping	Ice bag