

Anti-Padi4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Um07

Item Number: CM0021229



Introduction

High-quality anti-PADI4 rabbit polyclonal antibody validated for Western blot (WB), immunofluorescence (IF/ICC), and ELISA applications. Recognizes human, mouse, and rat PADI4. Suitable for studying citrullination, NETosis, and chromatin remodeling in immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PADI4 Rabbit pAb
Remarks/Alias	PAD; PAD4; PDI4; PDI5; PADI5; PADI4
Species	Human
GeneID (Human)	23569
GeneID	23569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-663 of human PADI4 (NP_036519.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UM07
Protein Weight	74kDa
Shipping	Ice bag