

Anti-Paep Rabbit Polyclonal Antibody For Wb, Elisa - P09466

Item Number: CM0030995



Introduction

CamelBio's Anti-PAEP rabbit polyclonal antibody detects human Glycodelin (P09466) via WB and ELISA. It supports research into fertilization, immunomodulation, and reproductive tissue studies.

[Learn More](#)

Parameter	Value
Product Name	PAEP Rabbit pAb
Remarks/Alias	GD; GdA; GdF; GdS; PEP; PAEG; PP14; ZIF-1; PAEP
Species	Human
GeneID (Human)	5047
GeneID	5047
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-180 of human PAEP (NP_002562.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P09466
Protein Weight	21kDa
Shipping	Ice bag