


Anti-Paf1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q8N7H5

Item Number: CM0009040



Introduction

Rabbit monoclonal antibody targeting human PAF1 (PD2), suitable for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for studying PAF1 complex function in transcription regulation and embryonic stem cell pluripotency.

[Learn More](#)

Parameter	Value
Product Name	PAF1 Rabbit mAb
Remarks/Alias	PD2; F23149_1; PAF1
Species	Human
GeneID (human)	54623
GeneID	54623
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PAF1 (NP_061961.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N7H5
Protein Weight	60kDa
Shipping	Ice bag