

Anti-Pah Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - P00439

Item Number: CM0002417



Introduction

Rabbit monoclonal antibody against human PAH (Phenylalanine-4-hydroxylase). Suitable for WB, IF/ICC, IF-P, ELISA. Cross-reacts with human, mouse, rat PAH. Ideal for phenylketonuria research and enzyme activity studies.

[Learn More](#)

Parameter	Value
Product Name	PAH Rabbit mAb
Remarks/Alias	PH; PKU; PKU1; PAH
Species	Human
Gene ID	5053
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PAH (P00439)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00439
Protein Weight	52kDa
Shipping	Ice bag