


## Anti-Pah Polyclonal Antibody For Wb, Elisa - P00439

Item Number: CM0022428



### Introduction

Rabbit polyclonal antibody targeting human PAH (phenylalanine-4-hydroxylase). Validated for WB and ELISA, with cross-reactivity to rat. Recombinant protein immunogen covering amino acids 110-452. Suitable for metabolism and genetic disorder research.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | PAH Rabbit pAb                                                                                        |
| Remarks/Alias    | PH; PKU; PKU1; PAH                                                                                    |
| Species          | Human                                                                                                 |
| GeneID (Human)   | 5053                                                                                                  |
| GeneID           | 5053                                                                                                  |
| Immunogen        | Recombinant protein corresponding to a sequence within amino acids 110-452 of human PAH (NP_000268.1) |
| Source           | Rabbit                                                                                                |
| Category         | Polyclonal Antibodies                                                                                 |
| Application      | WB, ELISA                                                                                             |
| Cross Reactivity | Rat                                                                                                   |
| SWISS            | P00439                                                                                                |
| Protein Weight   | 52kDa                                                                                                 |
| Shipping         | Ice bag                                                                                               |