

Anti-Pah Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P00439

Item Number: CM0016973



Introduction

Rabbit polyclonal antibody against phenylalanine-4-hydroxylase (PAH). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for phenylalanine metabolism, PKU research. Detects ~52kDa protein.

[Learn More](#)

Parameter	Value
Product Name	PAH Rabbit pAb
Remarks/Alias	PH; PKU; PKU1; PAH
Species	Human
GeneID (human)	5053
GeneID	5053
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PAH (NP_000268.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00439
Protein Weight	52kDa
Shipping	Ice bag