

Anti-Pai-1/Serpin E1 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P05121

Item Number: CM0027955



Introduction

Rabbit polyclonal antibody targeting human PAI-1/Serpin E1 (SERPINE1). Suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for studying fibrinolysis and cell adhesion.

[Learn More](#)

Specification	Details
Product Name	PAI-1/Serpin E1 Rabbit pAb
Remarks/Alias	PAI; PAI1; PAI-1; PLANH1; PAI-1/Serpin E1
Species	Human
Gene ID (Human)	5054
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 24-199 of human PAI-1/Serpin E1 (NP_000593.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05121
Protein Weight	45 kDa
Shipping	Ice bag