


Anti-Paics Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P22234

Item Number: CM0025307



Introduction

High-quality rabbit polyclonal antibody to human PAICS (bifunctional purine biosynthesis enzyme). Validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Cross-reacts with mouse and rat. Ideal for studying nucleotide metabolism and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	PAICS Rabbit pAb
Remarks/Alias	ADE2, AIRC, PAIS, ADE2H1, PAICSD, PAICS
Species	Human
Gene ID (Human)	10606
Gene ID	10606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 176-425 of human PAICS (NP_001072992.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22234
Protein Weight	47 kDa
Shipping	Ice bag