

Anti-Paip1 Polyclonal Antibody For Wb, Elisa - Q9H074

Item Number: CM0026007



Introduction

PAIP1 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, detecting human, mouse, and rat PAIP1. It aids research on translation initiation and mRNA turnover — Q9H074.

[Learn More](#)

Parameter	Value
Product Name	PAIP1 Rabbit pAb
Remarks/Alias	PAIP1
Species	Human
Gene ID	10605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PAIP1 (NP_877590.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H074
Protein Weight	54kDa
Shipping	Ice bag