


Anti-Pak1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13153

Item Number: CM0029234



Introduction

Rabbit polyclonal antibody against PAK1 (Q13153) validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for studying cytoskeleton dynamics, cell signaling, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	PAK1 Rabbit pAb
Remarks/Alias	IDDMSSD; p65-PAK; PAKalpha; alpha-PAK; PAK1
Species	Human
GeneID (human)	5058
GeneID	5058
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PAK1 (Q13153).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13153
Protein Weight	61kDa
Shipping	Ice bag