

Anti-Pak1p1 Polyclonal Antibody For Wb, Elisa - Q9Nwt1

Item Number: CM0021925

Introduction

Rabbit polyclonal antibody to PAK1IP1, a PAK1 kinase inhibitor and ribosome biogenesis factor. Validated for WB and ELISA with human and mouse cross-reactivity. Suitable for nucleolar signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PAK1IP1 Rabbit pAb
Remarks/Alias	PIP1; MAK11; WDR84; hPIP1; bA421M1.5; PAK1IP1
Species	Human
GeneID (Human)	55003
GeneID	55003
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-392 of human PAK1IP1 (NP_060376.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NWT1
Protein Weight	44kDa
Shipping	Ice bag