


Anti-Pak1/Pak2/Pak3 Rabbit Polyclonal Antibody For Wb, Elisa - Q13153/Q13177/O75914

Item Number: CM0020046



Introduction

Rabbit polyclonal antibody targeting PAK1, PAK2, and PAK3 (PAKalpha). Validated for WB and ELISA with reactivity in human, mouse, and rat tissues. Ideal for studying cytoskeleton dynamics and cell signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PAK1/PAK2/PAK3 Rabbit pAb
Remarks/Alias	PAKalpha; PAK1/PAK2/PAK3
Species	Human
Gene ID (Human)	5058/5062/5063
Gene ID	5058/5062/5063
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human PAK1/PAK2/PAK3 (NP_002567.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13153/Q13177/O75914
Protein Weight	65kDa
Shipping	Ice bag