

Anti-Pak2 Polyclonal Antibody For Wb And Elisa - Q13177

Item Number: CM0026202



Introduction

PAK2 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, and rat PAK2 (~58 kDa). Ideal for kinase signaling, apoptosis, and cytoskeleton research applications.

[Learn More](#)

Parameter	Value
Product Name	PAK2 Rabbit pAb
Remarks/Alias	KNO2; PAK65; PAKgamma; PAK2
Species	Human
GeneID (Human)	5062
GeneID	5062
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-240 of human PAK2 (NP_002568.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13177
Protein Weight	58kDa
Shipping	Ice bag