

## Anti-Pak2 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q13177

Item Number: CM0007867



### Introduction

High-quality rabbit monoclonal antibody against human PAK2 (Q13177), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for studying cytoskeleton regulation, apoptosis, and cell proliferation pathways.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | PAK2 Rabbit mAb                                                                                                    |
| Remarks/Alias    | KNO2; PAK65; PAKgamma; p58; PAK-2; S6/H4 kinase; PAK2                                                              |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 5062                                                                                                               |
| GeneID           | 5062                                                                                                               |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PAK2 (Q13177). |
| Source           | Rabbit                                                                                                             |
| Category         | Monoclonal Antibodies                                                                                              |
| Application      | WB, IF/ICC, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q13177                                                                                                             |
| Protein Weight   | 58kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |