


Anti-Pak3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75914

Item Number: CM0021955



Introduction

Rabbit polyclonal antibody targeting PAK3, a serine/threonine kinase involved in neuronal development and synaptic plasticity. Suitable for WB, IF/ICC, and ELISA applications. Recognizes PAK3 from human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	PAK3 Rabbit pAb
Remarks/Alias	ARA; bPAK; MRX30; MRX47; OPHN3; PAK-3; XLID30; PAK3beta; beta-PAK; PAK3
Species	Human
GeneID (Human)	5063
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75914
Protein Weight	62kDa
Shipping	Ice bag