

Anti-Pak4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O96013

Item Number: CM0015673



Introduction

Rabbit polyclonal antibody recognizing Human, Mouse PAK4 (64 kDa). Suitable for WB, IF/ICC, ELISA. Ideal for studying PAK4 in cytoskeleton regulation, cell migration, proliferation, and survival.

[Learn More](#)

Parameter	Value
Product Name	PAK4 Rabbit pAb
Remarks/Alias	PAK4
Species	Human
GeneID (Human)	10298
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 450-550 of human PAK4 (NP_005875.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O96013
Protein Weight	64kDa
Shipping	Ice bag